



# HACKATHON PROBLEM DESCRIPTION



Date: 20-21 Sep 2019

Prof. Arbind Kumar Gupta  
Subhamoy Chatterjee




Dayananda Sagar Institutions


# Anatomy of Heart



# Motivation: Ischemia and Fibrosis




# Cardiac Axis Planes




## Imaging Methods-




- Ultrasound (High frequency Sound waves)
- Computed tomography/ CT(X-ray)
- Magnetic Resonance Imaging/MRI (Magnetic field)
- MRI preferred over CT because of no radiation dosage



# Cross-sectional view of LV and Segmented Myocardium



Myocarditis with no LGE



Images on left side show Left Ventricle (LV) in cross-sectional view. The image at left-bottom shows scars in myocardium

Myocardium segmented from the above image and extracted.

# Problem 1

## Identification of Myocardium from MRI Short Axis (SAX) Cine Image



# Myocardium Segmentation

Some of the ideas are

- ❑ Template matching
- ❑ Hough transform for circle (disc)
- ❑ Convolutional Hough transform
- ❑ Scale Invariant convolutional Hough transform

Challenges are

- ❑ Shape may not be circular
- ❑ Size of myocardium may change




You will be provided with

- ❑ SAX DICOM images

You will need

- ❑ Pydicom [python package] for reading DICOM images

DICOM=Digital Imaging and COmmunications in Medicine








## Problem 2




How best to present scar tissues (Bull's Eye)



# Bull's Eye View



# How The Output May Look Like?



You will be provided with-

- LGE MR images and corresponding myocardium masks in 'png' format

You will need to-

- Stack up and align the slices
- Detect Fibrosis in each of them
- Create 17 segment display with %fibrosis printed scaled to  $(100/17)\% \sim 6\%$


## Problem 3

Unfold the myocardium volume and  
create a 3D curved MPR image from that




# Prototype



3D structure of Myocardium



Curved Plane Projection



# How The Output May Look Like?



You will be provided with-

- LGE MR images and corresponding myocardium masks in 'png' format

You will need to-

- Stack up and align the slices
- Create the 3D myocardium
- Create the curved plane projection for myocardium
- Display cross-section at different radial positions and slice planes